

A100-870

Aspect 100RT Fluoroscopic Procedure Table,
with Motorized Elevation, Trendelenburg
Tilt and Lateral Roll



- Carbon fiber tabletop with low x-ray attenuation
 - Tabletop: 24" x 80" (61 x 203 cm)
 - Optional Tabletop: 22" x 80" (56 x 203 cm) with facial cutout
 - 62" (157.5 cm) unobstructed imaging area
 - Motorized table movements with hand control
 - Elevation 30.4" – 42.2" (77 – 107 cm)
 - Trendelenburg Tilt $\pm 15^\circ$
 - Lateral Roll $\pm 20^\circ$
 - Auto return-to-level button and positional memory
 - Locking swivel casters
 - 450 lb. (204 kg) patient weight capacity
 - Table weight = 325 lbs. (147 kg)
 - 120 VAC operation (230 VAC optional)
 - Battery back-up operation
- Accessories included:**
- Tabletop Pad
 - Patient Restraint Straps
 - Facial pad, crescent shape
 - Disposable covers for facial pad, box of 50



Aspect Series Tables... C-arm Imaging At Its Best

Aspect Series multi-purpose tables are uniformly designed to provide these advantages...

- Small base pedestal and very large unencumbered imaging areas: 59" – 66" (149.9 – 167.6 cm), depending on model

- Long, slim-profile wheel base and high-strength steel lift column for...

✓ *Excellent stability and safety*

✓ *Very low table weigh for easy maneuverability and transport: 250 – 355 lbs. (113 – 161 kg), depending on model*

✓ *Long base provides protection from collisions that can damage carbon fiber tabletop*

✓ *450 lb. (204kg) "true" patient weight capacity; based on real patient positioning and over the entire range of tabletop positions*

- Sleek, easy-to-clean base and covers

- High-quality, carbon fiber tabletops and tabletop pads for excellent fluoroscopic imaging with minimal x-ray attenuation

- Multiple tabletop sizes and shapes to satisfy application-specific clinical requirements

- Highest quality components for exceptional durability and reliability. 3 year parts warranty.

- 120V AC operation (230V optional) with internal battery back-up operation. Less than 10 amps required.

- Certifications: ETL listed to UL 2601-1, CAN/CSA C22.2 No. 601.1-M90, EN 60601-1, IEC 60601-2-32.



(Table shown with standard rectangular tabletop)



(Table shown with optional tabletop with facial cutout)



Made in USA. Technical Specifications subject to change.



Product Specifications



Model Number:	100	100E	100R	100T	100RT	100RTL
Elevate (tabletop to floor)	Fixed	26 - 37.8" 66 - 96 cm	28.5 - 40.3" 72 - 102 cm	28 - 39.8" 71 - 101 cm	30.4 - 42.2" 77 - 107 cm	31.2 - 41.2" 79 - 104.6 cm
Trendelenburg Tilt	NA	NA	NA	±15°	±15°	±15°
Lateral Roll	NA	NA	20°	NA	20°	20°
Longitudinal Tabletop Motion	NA	NA	NA	NA	NA	10"
Standard Tabletop Size	24" x 80" 61 - 203 cm	24" x 80" 61 - 203 cm	24" x 80" 61 - 203 cm	24" x 80" 61 - 203 cm	24" x 80" 61 - 203 cm	24" x 80" 61 - 203 cm
Optional Tabletop Size	22" x 80" 56 - 203 cm	22" x 80" 56 - 203 cm	22" x 80" 56 - 203 cm	22" x 80" 56 - 203 cm	22" x 80" 56 - 203 cm	22" x 80" 56 - 203 cm
Imaging Area	66" 167.6 cm	66" 167.6 cm	62" 157.5 cm	62" 157.5 cm	62" 157.5 cm	59" 149.9 cm
Patient Capacity	450 lbs. 204 kg	450 lbs. 204 kg	450 lbs. 204 kg	450 lbs. 204 kg	450 lbs. 204 kg	450 lbs. 204 kg
Table Weight	250 lbs. 113 kg	290 lbs. 132 kg	325 lbs. 147 kg	325 lbs. 147 kg	325 lbs. 147 kg	355 lbs. 161 kg
Standard (●) and Optional (○) Accessories						
Tabletop Pad	●	●	●	●	●	●
Hand Control	NA	●	●	●	●	●
Foot Control	NA	○	○	○	○	○
ID410 Tabletop Width Extender Board	Optional: Available for 22" wide tabletop only.					
Clamp-on Accessory Rail	○	○	○	○	○	○
Armboard, rail mounted	○	○	○	○	○	○
Armboard, shoulder mount	○	○	○	○	○	○
Facial Support Pad	○	●	●	●	●	●
Restraint Straps	●	●	●	●	●	●
Table-side Radiation Shield	○	○	○	○	○	○

Tabletop Options



Standard Top

Cat. No. X100-995
24" x 80"
(61 x 203 cm)



Optional Top

Cat. No. X100-733
22" x 80"
(56 cm x 203 cm)
with facial cutout